## INFUSOR LV Portable Elastomeric Infusion System

| REF     | Description    | Nominal Vol. | Nominal Flow Rate* | Nominal Delivery Time | Minimum Vol. | Maximum Vol. | Approximate Residual Volume |
|---------|----------------|--------------|--------------------|-----------------------|--------------|--------------|-----------------------------|
| 2C2087K | INFUSOR LV 1.5 | 252 ml       | 1.5 ml/h           | 168 h                 | 227 ml       | 300 ml       | 3 ml                        |
| 2C2008K | INFUSOR LV 2   | 240 ml       | 2 ml/h             | 120 h                 | 216 ml       | 300 ml       | 3 ml                        |
| 2C2009K | INFUSOR LV 5   | 240 ml       | 5 ml/h             | 48 h                  | 216 ml       | 300 ml       | 3 ml                        |
| 2C2156K | INFUSOR LV 7   | 272 ml       | 7 ml/h             | 39 h                  | 245 ml       | 300 ml       | 3 ml                        |
| 2C2063K | INFUSOR LV 10  | 240 ml       | 10 ml/h            | 24 h                  | 216 ml       | 300 ml       | 3 ml                        |

\*The device is designed to operate at nominal flow rate using 5% Dextrose (D5W). Use of diluents other than D5W will affect the flow rate of the device (see Operating Conditions).

### Description

The INFUSOR LV system is an ambulatory, disposable device that uses an elastomeric reservoir to infuse medication. When filled, the INFUSOR LV system operates with a sustained internal pressure. Contents are delivered through a particulate matter filter and a flow restrictor.

The INFUSOR LV system is designed to provide a continuous flow of medication over the course of infusion at the nominal flow rate with an accuracy of  $\pm 10\%$  when filled with a volume ranging from **the minimum fill volume to the nominal volume plus the residual volume**.

The following factors may further affect the flow rate accuracy: (1) **Fill volume**, (2) **Temperature**, (3) **Viscosity of the medication**, (4) **The difference in height between the Fill Port and the Distal End Luer Lock/Flow Restrictor** and (5) **Catheter size** (refer to **Operating Conditions**).

The INFUSOR LV system Housing contains UV blocking capabilities effective to 380 nm which block UVB, UVC and most UVA rays.

This product is not manufactured using DEHP plasticizer.

This product is not manufactured with natural rubber latex.

### Indications for Use

The intended use of the INFUSOR LV system includes the slow, continuous intravenous, intra-arterial, subcutaneous, or epidural administration of medications. It may also include the slow, continuous infusion of medications directly into an intraoperative site or subcutaneously for postoperative pain management.

### Epidural Administration

Epidural administration of analgesics is limited to use with indwelling catheters specifically indicated for either short-term or long-term analgesic epidural drug delivery.

To prevent infusion of drugs that are not indicated for epidural use, do not use administration sets that incorporate injection sites.

It is strongly recommended that an INFUSOR LV system used for epidural drug delivery be clearly differentiated from an INFUSOR LV system used for any other route of administration.

### Warnings:

- It is the responsibility of the healthcare provider to ensure the use of an appropriate device when administering infusions of critical medications (such as Morphine) and life sustaining medications whose unplanned cessation, interruption or overdose would likely lead to serious injury or death.
- Do not fill with less than the minimum fill volume or more than the maximum fill volume (see product information table above).
- Assess the patient prior to use to ensure that the patient is appropriate for this device.
- It is the responsibility of the user to assure that the medication is prepared and administered in accordance with the drug manufacturer's package insert.
- Improper medication dosage can result in patient harm, including medication toxicity.
- Improper use of the device can result in patient harm, including medication toxicity.
- Epidural administration of drugs other than those indicated for epidural use could result in serious injury to the patient.**
- Single-use only. **Do not refill or sterilize.** Reuse or reprocessing of a single use device may lead to contamination and compromised device function or diminished structural integrity.
- Literature data suggests a possible association between the use of catheters in joint spaces during continuous intra-articular infusions and the development of chondrolysis.
- Check for solution clarity prior to administration. Do not use unless solution is clear.
- Do not use if either the Sterility Protector Cap or the Winged Luer Cap is not in place.
- Sterile Fluid Path. Use aseptic techniques for preparation and administration.
- Do not use if the package has been previously opened or the device appears damaged.
- Always closely supervise children and pets in the vicinity of the device during storage and use.
- The device is not designed to deliver whole blood nor cellular blood components.

### Cautions:

- The INFUSOR LV system must be filled in accordance with the procedures described under **Directions for Filling and Priming**.
- Store the unfilled device at room temperature and avoid exposure to extreme temperature.
- Allow the filled INFUSOR LV system to reach room temperature before infusing. Do not apply a heating source to bring the device to room temperature.
- The device is not designed for rapid infusion of medications.

### Information for Patients and Healthcare Provider:

**It is the responsibility of the healthcare provider to instruct the patient and caregivers of the following:**

- To achieve the desired medication temperature, ensure that the Distal End Luer Lock/Flow Restrictor is **secured to the skin**.
- Ensure that the Fill Port and the Distal End Luer Lock/Flow Restrictor are positioned at **approximately the same height** during the infusion. Use of a carrying bag or case may assist with positioning of the device.
- Check the INFUSOR LV system as directed by your healthcare provider. If there is any problem with the INFUSOR LV system, contact your healthcare provider.
- Observe the movement of the Elastomeric Reservoir along the Infusion Progress Lines over time to confirm that the medication is flowing.
- Do not disconnect the INFUSOR LV system before infusion is complete unless directed by your healthcare provider.
- When the INFUSOR LV system approaches end of the infusion, Internal Protrusions will begin to appear through the Elastomeric Reservoir. When all of the Internal Protrusions are visible, the infusion is essentially complete.
- Discard the device after a single use according to local regulations.
- Always closely supervise children and pets in the vicinity of the device during storage and use.
- Do not apply a heating source to bring the device to room temperature.

### Operating Conditions:

**Caution: Please note that a single or combination of the following factors may impact the flow rate:**

- Impact of Fill Volume**
  - The INFUSOR LV system is designed to operate at the nominal flow rate when it is filled with a volume ranging from **the minimum fill volume to the maximum fill volume** listed in the product information table above.
  - Do not fill with less than the minimum fill volume or more than the maximum fill volume (see product information table above).** Filling the device with less than the minimum fill volume will result in increases in flow rate and reductions in delivery time.

### 2. Impact of Temperature

- The INFUSOR LV system is designed to operate at the nominal flow rate when the medication is at a temperature of 31.1°C (88°F) for 2C2009K and 2C2063K and 33.3°C (92°F) for 2C2087K, 2C2008K and 2C2156K. To achieve the desired medication temperature, the Distal End Luer Lock/Flow Restrictor must be secured to the skin. Flow rate will decrease approximately 2.3% per 1°C (1.8°F) decrease in temperature and will increase approximately 2.3% per 1°C increase in temperature.

- Extended storage of the filled device in a refrigeration environment may result in a decrease in flow rate and an increase in delivery time.

### 3. Impact of Viscosity of the Medication

- The INFUSOR LV system is designed to operate at the nominal flow rate using 5% Dextrose (D5W) as the diluent. There will be an approximate 10% increase in the nominal flow rate if 0.9% Sodium Chloride (NS) is used.

### 4. Impact of Head Height

- The INFUSOR LV system is designed to operate at the nominal flow rate when the Fill Port and the Distal End Luer Lock/Flow Restrictor are positioned at **approximately the same height**. Flow rate will decrease approximately 0.5% for every one inch (2.54 cm) the Fill Port is positioned below the Distal End Luer Lock/Flow Restrictor and will increase approximately 0.5% for every one inch (2.54 cm) the Fill Port is positioned above the Distal End Luer Lock/Flow Restrictor.

### 5. Impact of Catheter Size

- The INFUSOR LV system is designed to operate at the nominal flow rate when used with a 22 gauge (3 French) catheter. Use of a catheter smaller than 22 gauge (3 French) will decrease the flow rate. Always refer to the directions for use provided by the catheter manufacturer.

- Note: As a result of intra-arterial pressure, the use of the device in intra-arterial access will lead to a reduction in flow rate.

### Directions for Use:

#### Mixing and Use Information

Refer to the drug manufacturer's package insert for drug reconstitution and storage procedures.

#### Directions for Filling and Priming

**The fluid path is sterile and nonpyrogenic if the Sterility Protector Cap and the Winged Luer Cap are in place. Do not use if the package has been previously opened or the device appears damaged.**

**Warning: Use aseptic technique throughout entire procedure.**

**Note: Failure to properly prime the device may result in delay of therapy.**

When withdrawing medication from a glass ampule, use a filter needle and remove the needle before filling the device.

**Note: Never use a needle when filling the INFUSOR LV system.**

- Remove the paper tape from the Tubing.
- Remove the Sterility Protector Cap from the INFUSOR LV system. Retain the Sterility Protector Cap to replace when filling is complete.
- Ensure that all air is removed from the syringe or filling device prior to filling. When possible, introduce the appropriate amount of diluent prior to filling with drug.
- Insert Luer tip of filled syringe or filling device into the Fill Port until secure. Do not over-tighten to avoid damaging the Fill Port.
- Fill the INFUSOR LV system.
- Fill the device with a volume ranging from the minimum fill volume to the maximum fill volume (see product information table above).
- Remove the syringe or filling device from the Fill Port.
- Replace and tighten the Sterility Protector Cap onto the Fill Port. Do not over-tighten to avoid damaging the Fill Port.
- A small air bubble may be visible within the Elastomeric Reservoir and is a normal occurrence. Do not attempt to remove the bubble.
- Prime the INFUSOR LV system by removing the Winged Luer Cap. Retain the Winged Luer Cap for future use.
- The INFUSOR LV system will prime slowly.
- Ensure that all air has been purged from the Tubing.
- If air bubbles are present in the Tubing, allow the device to flow until all air has been removed.
- If medication is not flowing at all to the Distal End Luer Lock/Flow Restrictor, force prime as follows: first, attach a Luer adapter (or stopcock) to the Distal End Luer Lock/Flow Restrictor. Second, attach a syringe (10 ml preferred) to the Luer adapter (or stopcock) and pull back the syringe plunger to apply continuous suction until fluid is observed in the syringe. Third, disconnect the syringe and Luer adapter (or stopcock) from the Distal End Luer Lock/Flow Restrictor.
- Visually confirm that the product is primed by observing two consecutive drops of medication from the tip of the Distal End Luer Lock/Flow Restrictor.
- Replace the Winged Luer Cap onto the Distal End Luer Lock/Flow Restrictor. Ensure that the Winged Luer Cap is securely connected after filling and priming. Over-tightening may cause the Winged Luer Cap to be difficult to remove.
- Label the INFUSOR LV system in accordance with applicable regulations and local institutional policy.
- Place the filled INFUSOR LV system in its pouch for transport to the patient.

### Directions for Administration:


**Refer to the Information for Patients and Healthcare Provider section.**

**Warning: Use aseptic technique throughout entire procedure.**

- To start the infusion, remove the Winged Luer Cap from the Distal End Luer Lock/Flow Restrictor. Visually confirm that medication is flowing by observing two consecutive drops of medication from the tip of the Distal End Luer Lock/Flow Restrictor before use.
- Ensure correct administration route. Securely connect the Distal End Luer Lock/Flow Restrictor to the catheter hub. Do not over-tighten.
- To achieve the desired medication temperature, ensure that the Distal End Luer Lock/Flow Restrictor is **secured to the skin**.
- Ensure that the Fill Port and the Distal End Luer Lock/Flow Restrictor are positioned at **approximately the same height** during infusion. Use of a carrying bag or case may assist with positioning of the device.

Baxter and Infusor are trademarks of Baxter International Inc.

## INFUSOR LV System Typical Flow Profile

(eng) **INFUSOR LV System Typical Flow Profile**  
 (fre) **Système INFUSOR LV - Profil de débit typique**  
 (ger) **Typisches Flussprofil des INFUSOR LV-Systems**  
 (spa) **Perfil característico del flujo del sistema INFUSOR LV**  
 (ita) **Tipico profilo di flusso del sistema INFUSOR LV**  
 (dut) **Typisch inlooppriofiel van het INFUSOR LV-systeem**  
 (swe) **INFUSOR LV-systemets typiska flödesprofil**  
 (nor) **Typisk infusjonsprofil for INFUSOR LV-systemet**  
 (dan) **Typisk flowprofil for INFUSOR LV systemet**  
 (fin) **INFUSOR LV -järjestelmän tyypillinen virtausprofiili**  
 (por) **Perfil de débito típico do sistema INFUSOR LV**  
 (gre) **Τυπικό προφίλ ροής συστήματος INFUSOR LV**  
 (chi) **INFUSOR LV 系統典型流速圖**

(eng) **Flow Rate**  
 (fre) **Débit**  
 (ger) **Flussrate**  
 (spa) **Velocidad de flujo**  
 (ita) **Velocità di flusso**  
 (dut) **Inloopsnelheid**  
 (swe) **Flödes hastighet**  
 (nor) **Infusjonshastighet**  
 (dan) **Flowhastighed**  
 (fin) **Virtausnopeus**  
 (por) **Taxa de débito**  
 (gre) **Ρυθμός ροής**  
 (chi) **流速**

(eng) **Delivery Time**  
 (fre) **Durée de l'administration**  
 (ger) **Verabreichungsdauer**  
 (spa) **Tiempo de administración**  
 (ita) **Tempo di somministrazione**  
 (dut) **Toedieningstijd**  
 (swe) **Infusionstid**  
 (nor) **Infusjonstid**  
 (dan) **Infusionstid**  
 (fin) **Infusion kesto**  
 (por) **Tempo de administração**  
 (gre) **Χρόνος χορήγησης**  
 (chi) **輸注時間**

(eng) Typical Flow Rate  
 (fre) Débit typique  
 (ger) Typische Flussrate  
 (spa) Velocidad de flujo establecida  
 (ita) Velocità di flusso tipica  
 (dut) Typische inloopsnelheid  
 (swe) Typisk flödes hastighet  
 (nor) Typisk infusjonshastighet  
 (dan) Typisk flowhastighed  
 (fin) Tyypillinen virtausnopeus  
 (por) Taxa de débito típica  
 (gre) Τυπικός ρυθμός ροής  
 (chi) 典型流速

(eng) Nominal Flow Rate  
 (fre) Débit nominal  
 (ger) Nominale Flussrate  
 (spa) Velocidad de flujo nominal  
 (ita) Velocità di flusso nominale  
 (dut) Nominale inloopsnelheid  
 (swe) Nominell flödes hastighet  
 (nor) Nominell infusjonshastighet  
 (dan) Nominel flowhastighed  
 (fin) Nimellinen virtausnopeus  
 (por) Taxa de débito nominal  
 (gre) Ονομαστικός ρυθμός ροής  
 (chi) 標定流速

(eng) ±10% of nominal flow rate over 100% of actual delivery time  
 (fre) ± 10 % du débit nominal au-delà des 100 % de la durée réelle de l'administration  
 (ger) ± 10 % der nominalen Flussrate während 100 % der tatsächlichen Verabreichungszeit  
 (spa) ±10 % de la velocidad de flujo nominal durante el 100 % del tiempo de administración real  
 (ita) ±10% della velocità di flusso nominale per il 100% del tempo di somministrazione effettivo  
 (dut) ±10% van de nominale inloopsnelheid bij 100% van de feitelijke toedieningstijd  
 (swe) ±10 % av den nominella flödes hastigheten över 100 % av den faktiska infusionstiden  
 (nor) ±10 % av nominell infusjonshastighet over 100 % av faktisk infusjonstid  
 (dan) ±10 % af den nominelle flowhastighed i løbet af 100 % af den faktiske infusionstid  
 (fin) ±10 % nimellisestä virtausnopeudesta yli 100 % todellisesta infuusion kestosta  
 (por) ±10% da taxa de débito nominal ao longo de 100% do tempo de administração real  
 (gre) ±10% του ονομαστικού ρυθμού ροής πάνω από το 100% του πραγματικού χρόνου χορήγησης  
 (chi) 在 100% 的實際輸注時間內，標定流速的 ±10%

| REF     | Description    | Nominal Vol. | Débit nominal* | Durée nominale de l'administration | Minimum Vol. | Maximum Vol. | Volume résiduel approximatif |
|---------|----------------|--------------|----------------|------------------------------------|--------------|--------------|------------------------------|
| 2C2087K | INFUSOR LV 1.5 | 252 ml       | 1,5 ml/h       | 168 h                              | 227 ml       | 300 ml       | 3 ml                         |
| 2C2008K | INFUSOR LV 2   | 240 ml       | 2 ml/h         | 120 h                              | 216 ml       | 300 ml       | 3 ml                         |
| 2C2009K | INFUSOR LV 5   | 240 ml       | 5 ml/h         | 48 h                               | 216 ml       | 300 ml       | 3 ml                         |
| 2C2156K | INFUSOR LV 7   | 272 ml       | 7 ml/h         | 39 h                               | 245 ml       | 300 ml       | 3 ml                         |
| 2C2063K | INFUSOR LV 10  | 240 ml       | 10 ml/h        | 24 h                               | 216 ml       | 300 ml       | 3 ml                         |

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**Baxter Healthcare Corporation**